

SURE-FLO INSUFFLATION TUBING

Insufflation Tubing
with .1 Micron ULPA Filter

10' of DEHP FREE Soft Tubing



2 Rotating Luer Locks
TM1TF10
10/Box



Tapered Connector
TM6000
10/Box



Quick Connector
TM102
10/Box

Insufflation Filter Sets
that allow for high volume
at low pressure settings.

LATEX FREE & DEHP FREE

Filters Particulate Matter
of .1 Micron With
99.997% Efficiency

Hydrophobic ULPA Filtration

15/22 mm Connector
TM1522
10/Box



SUCTION/IRRIGATION



Irrigation Pump

- No Batteries
- 1 Liter and Dual 1 Liter Capacity
- Noise Free
- Locking Pressure Adjustment
- Fast Loading & Unloading

"Green Initiative"
'Eliminate Use of Batteries'



Suction Irrigation Set

- No Batteries
- Dual Spike
- Tubing Split (Irrigation/Suction)
- Multiple Probes & Cautery Probes
- Ergonomic Hand Piece



P.O. Box 800192
Santa Clarita, Ca 91380-0192
1.800.971.9066 or 661.257.0788
Fax: 661.257.0688
www.bhsurgical.com